
Liquid Cell Holder

The X-Rite Liquid Cell Holder is designed to hold a Starna Type 96 rectangular cell for liquid measurement. The holder features a black chamber with a white ceramic backing that encloses the liquid-filled cell. A cell lid is provided to prevent spilling. With the Liquid Cell Holder, we recommend holding the instrument in a vertical position using the X-Rite Benchtop Stand (sold separately).

For use with	Part number
ALL SPs and 900 series	DRS-810
Replacement Cell Lid	SP68-96G10L

Cell Holder for Powder or Liquid

The X-Rite Cell Holder is designed to hold a Hellma Circular Cell (Hellma part number 692.091) for liquid or powder measurement. The fixture attaches to the shoe of an X-Rite portable spectrophotometer and includes a Hellma Circular Cell, which is enclosed in a black chamber to prevent ambient light from affecting the measurement. With the Liquid Cell Holder, we recommend holding the instrument in an upside-down position using the X-Rite Benchtop Stand (sold separately).

For use with	Part number
SP60, SP62 and SP64	DRS-811
SP68 though SP98	SP68-811
900 Series	968-811

Universal Sample Holder

